

Product Contents

Alkaline Phosphatase, Calf Intestinal (CIAP):

Part No.	Conc. (μM)	Size (units)
M182A	1	1,000
M282A	20	1,000

Part# 9PIM182
Revised 9/98

Description: Alkaline Phosphatase catalyzes the hydrolysis of 5'-phosphate groups from DNA, RNA and both ribo- and deoxyribonucleoside triphosphates.

Enzyme Storage Buffer: Alkaline Phosphatase, Calf Intestinal (CIAP), is supplied in 10mM Tris-HCl (pH 8.0), 1mM MgCl_2 , 0.1mM ZnCl_2 , 50mM KCl and 50% glycerol.

Alkaline Phosphatase 10X Reaction Buffer (M183A): When the 10X Reaction Buffer supplied with this enzyme is diluted 1:10, it has a composition of 50mM Tris-HCl (pH 9.3 at 25°C), 1mM MgCl_2 , 0.1mM ZnCl_2 and 1mM spermidine.

Source: Calf intestinal mucosa.

Unit Definition: One unit is defined as the amount of enzyme required to catalyze the hydrolysis of 1 μmol of 4-nitrophenyl phosphate per minute at 37°C in 1M diethanolamine, 10.9mM 4-nitrophenyl phosphate, 0.5mM MgCl_2 (pH 9.8). See the unit concentration on the Product Information Label.

Storage Temperature: For long term storage (infrequent use; 1-2 times per month), store at -70°C. For daily/weekly use, store at -20°C. Avoid multiple freeze-thaw cycles. See the expiration date on the Product Information Label.

Quality Control Assays

Contaminant Activity

Endonuclease Assay: To test for endonuclease activity, 1 μg of Type I supercoiled plasmid DNA is incubated with 5 units of Calf Intestinal Alkaline Phosphatase in 1X Reaction Buffer for one hour at 37°C. Following incubation, the supercoiled DNA is visualized on an ethidium bromide-stained agarose gel to verify the absence of visible nicking or cutting.

DNase and RNase Assay: To test for nuclease activity, 50ng of radiolabeled DNA or radiolabeled RNA is incubated with 5 units of Calf Intestinal Alkaline Phosphatase in 1X Reaction Buffer for one hour at 37°C, and the release of radiolabeled nucleotides is monitored by scintillation counting of TCA-soluble material. Minimum passing specification is $\leq 3\%$ release for DNase and $\leq 3\%$ release for RNase.

Blue/White Assay: Five micrograms of pGEM[®]-3Zf(+) Vector is digested with representative restriction enzymes (leaving 5'-termini, 3'-termini or blunt ends). The termini are treated with 5 units of Calf Intestinal Alkaline Phosphatase for 1 hour at 37°C, kinased and ligated. The religated plasmid is then transformed into JM109 cells which are plated on X-Gal/PTG/Amp₁₀₀ plates. A minimum of 200-400 colonies are counted. White colonies result from transformation with ligated plasmids with damaged ends. These white colonies represent the number of false positives expected in a typical cloning experiment. Enzymes which generate overhangs, such as *EcoR* I, *Hind* III or *Kpn* I, must produce fewer than 2% white colonies, and blunt cutting enzymes, such as *Hinc* II, must produce fewer than 5% white colonies (1). Transformation efficiency must be $\geq 1 \times 10^5$ cfu/ μg DNA.

Reference

- Hung, L. *et al.* (1991) A blue/white cloning assay for quality control of DNA restriction and modifying enzymes. *Promega Notes* 33, 12.



Promega

Promega Corporation

2800 Woods Hollow Road	
Madison, WI 53711-5399	USA
Telephone	608-274-4330
Toll Free	800-356-9526
Fax	608-277-2516
Internet	www.promega.com
Europe	www.euro.promega.com

PRODUCT USE LIMITATIONS, WARRANTY, DISCLAIMER

Promega manufactures products for a number of intended uses. Please refer to the product label for the intended use statements for specific products. Promega products contain chemicals which may be harmful if misused. Due care should be exercised with all Promega products to prevent direct human contact.

Each Promega product is shipped with documentation stating specifications and other technical information. Promega products are warranted to meet or exceed the stated specifications. Promega's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Promega makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, PRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO PROMEGA PRODUCTS. In no event shall Promega be liable for claims for any other damages, whether direct, incidental, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of Promega products to perform in accordance with the stated specifications.

*U.S. Pat. No. 4,766,072 has been issued to Promega Corporation for transcription vectors having two different bacteriophage RNA polymerase promoter sequences separated by a series of unique restriction sites into which foreign DNA can be inserted.

© 1997, 1998 Promega Corporation. All Rights Reserved.

pGEM is a trademark of Promega Corporation and is registered with the U.S. Patent and Trademark Office.

All specifications are subject to change without prior notice.

Product claims are subject to change. Please contact Promega Technical Services or access the Promega online catalog for the most up-to-date information on Promega products.

Part# 9PIM182
Printed in USA. Revised 9/98

I. Description

Calf intestinal alkaline phosphatase (CIAP) catalyzes the hydrolysis of 5'-phosphate groups from DNA, RNA, and ribo- and deoxyribonucleoside triphosphates. This enzyme is used to prevent recircularization and religation of linearized cloning vehicle DNA by removing phosphate groups from both 5'-termini (1-5).

II. Reaction Conditions

A. Dephosphorylation of 5' Overhangs

Reagents to Be Supplied by the User

(Solution compositions are provided in Section III.)

- 10mM Tris-HCl (pH 8.0)
- CIAP stop buffer
- TE-saturated phenol:chloroform
- chloroform:isoamyl alcohol (24:1)
- 7.5M ammonium acetate (pH 5.5)
- ethanol, 100% and 70%

1. Dilute sufficient CIAP for immediate use in CIAP 1X Reaction Buffer to a final concentration of 0.01u/μl. Each pmol of DNA ends will require 0.01u CIAP. (1μg of 1,000bp DNA = 1.52pmol DNA = 3.03pmol of ends.)
2. Purify the DNA to be dephosphorylated by ethanol precipitation and resuspend the pellet in 40μl of 10mM Tris-HCl (pH 8.0). Set up the following reaction:

DNA (up to 10 pmol of 5'-ends)	40μl
CIAP 10X Reaction Buffer	5μl
Diluted CIAP (0.01u/μl)	<u>up to 5μl</u>
	50μl

3. Incubate at 37°C for 30 minutes.
4. Add another aliquot of diluted CIAP (equivalent to the amount used in Step 2) and continue incubation at 37°C for an additional 30 minutes.
5. Add 300μl of CIAP stop buffer. Phenol:chloroform extract and ethanol precipitate by adding 0.5 volume 7.5M ammonium acetate (pH 5.5) and 2 volumes of 100% ethanol to the final aqueous phase.

Note: CIAP may be added directly to digested DNA. Add 5μl CIAP 10X Reaction Buffer, 0.01u CIAP/pmol of ends and deionized water to a final volume of 50μl (6).

B. Dephosphorylation of 5' Recessed or Blunt Ends

When 5' recessed or blunt end DNA fragments are used as substrate, incubate at 37°C for 15 minutes and then at 56°C for 15 minutes. Then add a second aliquot of CIAP and repeat the incubations at both temperatures. The higher temperature ensures accessibility of the recessed end (7).

III. Composition of Buffers and Solutions

CIAP stop buffer

10mM	Tris-HCl (pH 7.5)
1mM	EDTA (pH 7.5)
200mM	NaCl
0.5%	SDS

TE buffer

10mM	Tris-HCl (pH 8.0)
1mM	EDTA

TE-saturated phenol:chloroform

Mix equal parts of TE buffer and phenol and allow the phases to separate. Then mix 1 part of the lower, phenol phase with 1 part of chloroform:isoamyl alcohol (24:1).

IV. References

1. Sambrook, J., Fritsch, E.F. and Maniatis, T. (1989) *Molecular Cloning: A Laboratory Manual*, Cold Spring Harbor Laboratory, Cold Spring Harbor, New York.
2. Seeburg, P.H. *et al.* (1977) Nucleotide sequence and amplification in bacteria of structural gene for rat growth hormone. *Nature* **220**, 486.
3. Ullrich, A. *et al.* (1977) Rat insulin genes: construction of plasmids containing the coding sequences. *Science* **196**, 1313.
4. Meyerowitz, E.M. *et al.* (1980) A new high-capacity cosmid vector and its use. *Gene* **11**, 271.
5. Grosfeld, F.G. *et al.* (1981) Isolation of beta-globulin-related genes from a human cosmid library. *Gene* **13**, 227.
6. *Protocols and Applications Guide*, Third Edition (1996) Promega Corporation.
7. Perbal, B. (1988) In: *A Practical Guide to Molecular Cloning*, Second Edition, John Wiley and Sons, New York.